

NAME: \_\_\_\_\_

PHONE: \_\_\_\_\_

COMPANY/ORGANIZATION: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

DEALER/REP: \_\_\_\_\_

SALES REP: \_\_\_\_\_ QUOTE: \_\_\_\_\_



## CLASSIC SERIES MAJESTIC DELUXE ELECTRIC TABLE

### STANDARD FEATURES:

- Lifting Capacity: 600 lbs.
- Height range: 27.5-39.5"
- Dual action face rest base
- Salon head support and neck pillow
- Standard single foot pedal
- 9" on center dual-end accessory ports
- Holes for swivel side arm supports - head end
- Swivel side arm supports
- 3" Foam system and deluxe cushion
- Rounded corners
- Precision cut steel base with column lift pedestal
- Whisper quiet lift system
- Manually adjusting back section and knee sections



Sleek and modern, our Deluxe model in the Majestic Electric Massage Table line provides effortless yet incredibly comfortable support in every position. You can smoothly transition clients from supine to seated with the standard manually tilting back and knee sections. Its column base allows for generous leg and knee room, 12-inch height range adjustment and ease of operation, with just a tap of the toe. The Deluxe model includes a dual action face rest, salon head support and neck pillow, swivel side arms and rounded corners. Crescent pillow not included.

### STANDARD TABLES:

**FLO-2807** - Basic 28"  
28"W x 73"L x 27.5-39.5"H

Qty \_\_\_\_\_

**FLO-3007** - Basic 30"  
30"W x 73"L x 27.5-39.5"H

Qty \_\_\_\_\_

**FLO-3307** - Basic 33"  
33"W x 73"L x 27.5-39.5"H

Qty \_\_\_\_\_

### STANDARD UPHOLSTERY CHOICES:

- |   |   |
|---|---|
| <input type="checkbox"/>  Agate      | <input type="checkbox"/>  Hunter Green |
| <input type="checkbox"/>  Black      | <input type="checkbox"/>  Mushroom     |
| <input type="checkbox"/>  Burgundy   | <input type="checkbox"/>  Navy         |
| <input type="checkbox"/>  Buff       | <input type="checkbox"/>  Purple       |
| <input type="checkbox"/>  Chocolate  | <input type="checkbox"/>  Teal         |
| <input type="checkbox"/>  Flint Gray |   |

**OPTIONS:**

- NTU** - Natural Touch Upholstery
- NT/TRU** - Natural Touch Upholstery & Touch Response Memory Foam System
- BR** - Breast Recesses
- UP DPOS-B** - Upgrade to Deluxe Adjustable Face Rest Base
- UP PPOS-B** - Upgrade to Pivot Posi-tilt Face Rest Base
- Sub Hand** - Substitute Hand Control for Foot Control
- Add Foot** - Add Second Foot Control
- Add Hand** - Add Hand Control
- PTB** - Power Tilt Back with Hand Control
- Sub Foot PTB** - Substitute Foot Control for Hand Control
- Add Foot PTB** - Add Foot Control to Power Tilt Back
- Add Hand PTB** - Add Hand Control to Power Tilt Back
- CR-Pill** - Classic Series Crescent Pillow with Standard Upholstery
- CR-Pill NT** - Classic Series Crescent Pillow with Natural Touch Upholstery
- CR-Pill NT/TRU** - Classic Series Crescent Pillow with Natural Touch Upholstery and Touch Response Foam
- 8 OC** - Custom Configuration 8" on Center Face Rest Holes

**UPGRADED NATURAL TOUCH UPHOLSTERY CHOICES:**

- |   |   |
|---|---|
| <input type="checkbox"/>  Amethyst | <input type="checkbox"/>  Stainless Steel Gray |
| <input type="checkbox"/>  Sapphire | <input type="checkbox"/>  Sand Dune            |
| <input type="checkbox"/>  Jade     | <input type="checkbox"/>  White Tie            |
| <input type="checkbox"/>  Ivory    | <input type="checkbox"/>  Morning Dew          |
| <input type="checkbox"/>  Ebony    | <input type="checkbox"/>  Latte                |



*Breast Recesses Option - Breast/Scaupula Recesses include dual cutouts with central sternum/spine support (includes custom plugs)*



*PTB - Power Tilt Back with Hand Control*



*UP DPOS-B - Upgrade from Dual Action Face Rest Base to Deluxe Adjustable. Full range of adjustment is 4-6".*